

## Walking Aids

### Compact Indoor Rollator

This lightweight sturdy compact rollator with a height adjustable handle is perfect for indoor and around the home use. With a delta handle and brake it suits those with limited upper limb mobility, stroke survivors or limb amputees. The removable tray and storage bag makes it easy to transport meals comfortably, great for carrying utensils or personal belongings around the house. This rollator is easy to fold, ideal for storage and transporting.

- \* Lightweight sturdy frame
- \* Delta handle and brake
- \* Perfect for indoor use around the home
- \* Handy tray for carrying meals
- \* Removable storage bag
- \* 8" Castors
- \* Easy folding for storage and transport
- \* Adjustable handle height
- \* Weight capacity 136kg



Product Code: RM926G

#### Specifications

Tray size: 32x23x2.5 (cm)

Bag size: 39x27x14 (cm)

Overall Height	Overall Length	Overall Width	Folded Width	Unit Weight	Weight Capacity
85.5-96.5	75	58.5	42	6.7	136

Measurements in cm, Weight in kg