

Wheelchair

Compact Lite Wheelchair

The Compact Lite is as the name suggests compact all round great wheelchair. Foldable, light in weight makes it a great chair not just for the occupant, but also the carer. A smart, sporty and classy looking compact unit, it is comfortable and perfect for every day use.

- * Lightweight frame, easy to manage
- * Folding backrest
- * Attendant Brake
- * Swing up leg rests with 2 height options
- * AP - Solid 6" front castor, 12" rear wheel
- * SP - Solid 6" front castor, 20" rear wheel
- * Rear wheel locking brake
- * Seat belt included
- * Weight capacity 110kg



Product Code: RM836-AP



Product Code: RM836-SP

Specifications

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Rear Wheel	Unit Weight	Weight Capacity
RM836-SP	92	99	60	38	45	46	21	44	20"	12.2	110
RM836-AP	92	97	60	38	45	46	21	44	12"	10.9	110

Measurements in cm, Weight in kg

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com