

## Walking Aids *Bariatric*

### Deluxe Bariatric Rollator

The aluminium heavy duty deluxe rollator is suitable for indoor or outdoor use. It includes height adjustable handles and height adjustable frame providing a seat option to suit shorter to taller stature. A soft comfortable seat with a detachable shopping bag underneath. It has 8" PVC wheels and a weight capacity of 181kgs.

- \* Heavy duty aluminium frame
- \* Oval tubular frame construction
- \* Height adjustable handles and frame
- \* Soft touch brake
- \* 8" PVS wheels
- \* Ergonomic hand grips for comfort
- \* Fabric padded seat with bag
- \* Weight user: 181kg



Product Code: VP187D

#### Specifications

| Overall Height | Overall Length | Overall Width | Seat Depth | Seat Width | Seat Height | Folded Width | Unit Weight | Weight Capacity |
|----------------|----------------|---------------|------------|------------|-------------|--------------|-------------|-----------------|
| 83-103         | 80             | 66            | 40         | 47         | 50-58       | 45           | 13          | 181             |

Measurements in cm, Weight in kg