

## Bedroom

### Deluxe Overbed Table

The deluxe Multi Purpose Table is a convenient, mobile table for use around the home, particularly useful as an overbed table. The main table top able can be set to different angles from flat to 45 deg. with an additional side table providing excellent space for a drink or other smaller items.



Product Code: VG832A

- \* Suitable for use over the bed or beside a chair
- \* Powder coated frame
- \* Beech melamine top with upstands
- \* Height adjustable
- \* Angled top can tilt from flat to approx. 45deg.
- \* Fixed side table for drinks
- \* 4 easy roll locking castors
- \* Maximum weight limit 20kg
- \* Clearance required 95mm

#### Specifications

Overall Height	Overall Length	Overall Width	Unit Weight
76-117	81	39	10

Measurements in cm, Weight in kg