

## Walking Aids

### Duo Transporter Rollator

The Duo Transporter is a high quality mobility aid having a dual purpose, converting from a rollator to a wheelchair in seconds. This allows the user to walk with support of the rollator and to sit down and rest when needed. The adjustable, removable foot rests allow it to be used as a wheelchair with a flip able backrest to enable the user to be pushed. Suitable for indoor and outdoor use it is provided with a motion and park brake, and height adjustable handle bars. The strong aluminium frame is easily foldable and great for frequent travellers.



- \* Aluminium frame
- \* Double cross brace frame
- \* Great for frequent travellers
- \* Height adjustable handle bars
- \* Motion and park brake
- \* Flip able padded backrest
- \* Easy to fold and transport
- \* Removable bag
- \* Max weight user: 110kg
- \* Colour: Gold

Product Code: RM199GD

#### Specifications

| Overall Height | Overall Length                | Overall Width | Seat Depth | Seat Width | Seat Height | Folded Width | Unit Weight | Weight Capacity |
|----------------|-------------------------------|---------------|------------|------------|-------------|--------------|-------------|-----------------|
| 85-97          | Rollator:74<br>Wheelchair:101 | 63            | 30         | 48         | 55.5        | 28           | 9.7         | 110             |

Measurements in cm, Weight in kg