

## Walking Aids

### Forearm Walker

The aluminium lightweight forearm rollator allows the user to support their weight on the arms cuffs to reduce strains on the wrists and hands. The forearm cuffs are contoured and height adjustable. Featuring 10" front and 8" rear PVC wheels, and includes a detachable shopping bag. Easily foldable and compact for storage.



Product Code: VP187F

- \* Aluminium lightweight frame
- \* Height adjustable handles
- \* Easily foldable and compact
- \* Includes fabric seat and detachable shopping bag
- \* 10" front and 8" rear PVC wheels
- \* Weight capacity 136kg

#### Specifications

Overall Height	Overall Length	Overall Width	Armrest Height	Seat Depth	Seat Width	Seat Height	Folded Width	Unit Weight	Weight Capacity
108-123	67	62	84-109	23	44	55.5	29.5	12	136

Measurements in cm, Weight in kg