

Wheelchair

G3 Wheelchair

The G3 Wheelchair is equipped with a number of special features along with a host of accessory choices, makes this a very desirable chair for everyday use. The self propel model is available with a standard brake or a drum brake along with care assist handbrake.



- * Sturdy aluminium frame
- * Angle and height adjustable seat
- * Folding backrest
- * Care assist brake
- * Quick release 24" mag rear wheels with PU Tyres
- * Swing away removable leg rests
- * Padded nylon upholstery
- * 7" front castor
- * Optional drum brake (S/P only)
- * Anti-tip wheels
- * Weight capacity 150kg

Optional Stump Support



Optional Elevating Leg Rest



Product Code:

RM138-46 S/P 46cm Seat

RM138-46D S/P 46cm Seat with Drum Brake

RM138-51 S/P 51cm Seat

RM138-51D S/P 51cm Seat with Drum Brake

Specifications

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Footrest Height	Unit Weight	Weight Capacity
RM138-46	74.5-90.5	107	67	41	46	46-50	22	38	17	16.9	150
RM138-51	74.5-90.5	107	71	41	51	46-50	22	38	17	17.1	150

Measurements in cm, Weight in kg