

## Wheelchair

### G3 Wheelchair

The G3 Wheelchair is equipped with a number of special features along with a host of accessory choices, makes this a very desirable chair for everyday use. The self propel model is available with a standard brake or a drum brake along with care assist handbrake.



- \* Sturdy aluminium frame
- \* Height adjustable seat
- \* Folding backrest
- \* Care assist brake
- \* Quick release 24" rear wheels with PU Tyres
- \* Swing away removable leg rests
- \* Padded nylon upholstery
- \* 7" front castor
- \* Optional drum brake (S/P only)
- \* Anti-tip wheels
- \* Weight capacity 150kg

Optional Stump Support



Optional Elevating Leg Rest



#### Product Code:

RM138-46 S/P 46cm Seat

RM138-46D S/P 46cm Seat with Drum Brake

RM138-51 S/P 51cm Seat

RM138-51D S/P 51cm Seat with Drum Brake

#### Specifications

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Footrest Height	Unit Weight	Weight Capacity
RM138-46	74.5-90.5	107	67	41	46	46-50	22	38	17	16.9	150
RM138-51	74.5-90.5	107	71	41	51	46-50	22	38	17	17.1	150

Measurements in cm, Weight in kg

## Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email [support@ausnewhomecare.com](mailto:support@ausnewhomecare.com)