

Wheelchair

G4 Plus Wheelchair

The G4 Wheelchair is multi-adjustable allowing the chair to be tailored to the individual. Equipped with a number of important features along with a host of accessory choices, makes this a very desirable chair for everyday use. This model has an optional drum brake providing extra comfort for the passenger.

- * Sturdy aluminium frame
- * Height adjustable seat
- * Height adjustable attendant assist brake
- * Height and angle adjustable backrest
- * Flip up, removable, height adjust, sliding, padded armrest
- * Swing away removable leg rests
- * Padded nylon upholstery
- * Quick release 24" mag rear wheels with PU Tyres
- * 7" front caster
- * Optional drum brake
- * Width extendable foot plate
- * Height adjustable Anti-tip bars
- * Weight capacity 140kg



Drum brake pictured

Product Code:

- RM144-46 (AWC145)** 46cm
- RM144-46D** 46cm with Drum Brake
- RM144-51 (AWC144)** 51cm
- RM144-51D** 51cm with Drum Brake



Optional Elevating Leg Rest



Optional Stump Support



Specifications

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Footrest Height	Unit Weight	Weight Capacity
RM144-46	100-107	106	67	40-50	46	51	21-28	48-53	17	17.4	140
RM144-51	100-107	106	72	40-50	51	51	21-28	48-53	17	17.4	140

Measurements in cm, Weight in kg