

Wheelchair

G6 Excel Bariatric Wheelchair

The G6 Excel Bariatric Wheelchair is multi-adjustable with individual tailored options allowing this chair to be modified to the individual user. A number of important features along with a host of accessory choices, makes this a very desirable chair for everyday use with the comfort of a padded seat. Care assist brake gives a little extra assurance for the career. This model is available with a standard brake or with a drum brake providing extra comfort for the passenger.



- * Strong steel frame
- * Angle and height adjustable seat
- * Swing away removable leg rests
- * Flip back and removable armrests
- * Quick release 24" rear wheels
- * 7" front caster
- * Attendant assist hand brake
- * Optional drum brake
- * Reclining backrest
- * Weight capacity 250kg



Optional 1 Piece
Footplate

Product Code:

RM160-56 56cm

RM160-56D 56cm with Drum Brake

RM160-61 61cm

RM160-61D 61cm with Drum Brake



Optional Stump
Support



Optional Elevating
Leg Rest

Specifications

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Footrest Height	Unit Weight	Weight Capacity
RM160-56	103-116	102	78	46	56	53	23	45	6-24	26.3	250
RM160-61	103-116	102	82	46	61	53	23	45	6-24	26.9	250

Measurements in cm, Weight in kg