

Wheelchair

G6 Excel Bariatric Wheelchair

The G6 Excel Bariatric Wheelchair is multi-adjustable with individual tailored options allowing this chair to be modified to the individual user. A number of important features along with a host of accessory choices, makes this a very desirable chair for everyday use with the comfort of a padded seat. Care assist brake gives a little extra assurance for the career. This model is available with a standard brake or with a drum brake providing extra comfort for the passenger.



- * Strong steel frame
- * Angle and height adjustable seat
- * Swing away removable leg rests
- * Flip back and removable armrests
- * Quick release 24" rear wheels
- * 7" front caster
- * Attendant assist hand brake
- * Optional drum brake
- * Reclining backrest
- * Weight capacity 250kg



Specifications 56/61

Overall height	103-116	Seat height.....	53
Overall depth	102	Backrest height	45
Overall width	78/82	Armrest height.....	23
Seat depth.....	46	Unit weight	26/27
Seat width.....	56/61	Weight capacity	250

Measurements in cm, weight

Product Code: RM160-56 56cm Seat
Product Code: RM160-56D Drum Brake
Product Code: RM160-61 61cm Seat
Product Code: RM160-61D Drum Brake

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com