

Walking Aids

Aidapt Knee Walker

The Innovative Knee Walker is a great option to help with the recovery from a foot injury, break or sprain, allowing the user to evenly distribute their body weight on the device. With the supplied cushioned knee pad to support injured leg, ankle or foot providing maximum comfort. It has a sturdy dual frame which adds stability and durability helping with easy steering and stopping.

- * Sturdy dual frame
- * Compact for easy transportation
- * Height adjustable handle bar and knee pad
- * Detachable basket
- * Fits adults, young adults and elderly.
- * Weight capacity 113kg



Product Code: VP186

Additional Sheepskin Knee Pad Cover available

Product Code: SS6000



Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Unit Weight	Weight Capacity
86-106	80	43	8	70	48-58	10	113

Measurements in cm, Weight in kg