



Fact Sheet

Bathroom



Medina Bath Board

The Medina Bath Board is a high quality, reinforced plastic bath board.

- * Modern design
- * Slatted surface provides adequate drainage
- * Wipe clean surface
- * Simple to install with self locking cleats
- * Both available with or without handle
- * Weight Capacity 159kg



Product Code: VR116A (L:600mm)
VR117A3 (L: 685mm)
VR118A (L:710mm)

Specifications

- Overall Depth - 4 slat (mm) 250
- Overall Depth—6 slat (mm) 350
- Seat Height (mm) 140
- Seat Width (mm) 685
- Width Adjustment 535-625
- Net weight 3.2kg



Product Code: VR117H



Product Code: VR125A (L:685mm)

Medina Bath Seat—Bariatric

The Medina Bath Seat provides extra height enabling easy access in and out of the bath.

- * Constructed from high quality, fully reinforced slatted plastic
- * Available in 2 heights (6" & 8")
- * Suitable for majority of bath's
- * Slatted for easy drainage
- * Modern easy to clean surface
- * Secured to bath with 4 large suction caps
- * Easily removed by pulling on the integral tabs
- * Finished in soft grey
- * Weight Capacity 190kg



Specifications

- Dimensions (6") D290xW450xH150
- Dimensions (8") D290xW450xH200
- Net weight 2.4kg

Product Code: VR127A (6")
VR128A (8")

Your Local Agent

QLD 1/14-16 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW 5/143 Industrial Rd, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au