

## Bathroom *Bariatric*

### Over Toilet Aid










The Bariatric 406 Over Toiler Aid constructed from aluminium oval tubing offers extra strength and security in its frame and support. The high density wider pvc seat offers ultimate comfort and comes with easily fitting modesty cover and splash guard. Fitted with large rubber suction tips to reduce slipping.

- \* Lightweight toughened powder coated aluminium frame
- \* Large suction tips to reduce slipping
- \* Wider seat for increased comfort
- \* 65cm between arms
- \* Modesty cover
- \* Splash guard
- \* Height adjustable
- \* Weight capacity 305kg



Product Code: RM406-65

### Specifications

overall height	overall depth	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
								
67-79	54	87	45	65	43-55	24	8.65	305

Measurements in cm, Weight in kg