

Walking Aids

Aidapt Rollator 8"

The Aidapt 8" Rollator is an affordable option, made up of a strong steel frame. With a comfortable padded seat, backrest and 8" solid castors, this rollator is suitable for daily use both indoors and out.

- * Strong steel frame
- * Height adjustable handles
- * Ergonomic handles for comfort
- * Padded seat and back
- * 8" solid castors
- * Max weight user: 130kg



Specifications

Overall height	790-900	Seat height	52
Overall depth	69	Width between arms	46
Overall width	61	Unit weight	8
Seat depth.....	32	Weight capacity	130
Seat width.....	34		

Measurements in cm, weight

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com