

Walking Aids

RM X Frame Rollator

The Deluxe X lightweight, aluminium four wheeled Rollator has a very modern and stylish design with many great features. Featuring a full flat fold design, makes it ideal not only for storage but also transportation. With anatomically-designed handgrips, arthritic friendly loop cable brakes and height adjustable handles, the Aidapt Deluxe X is suitable for the elderly, infirm or less able users with little grip power.











- * Lightweight Aluminium frame
- * Easy side fold
- * Height adjustable
- * Canvas seat with padded backrest
- * Walking stick holder
- * Detachable shopping bag
- * Max weight user: 130kg

Available in Silver



Product Code: RM200S

Specifications

overall height	overall depth	overall width	seat depth	seat width	seat height	folded depth	castor size	unit weight	weight capacity
									
78-90	69	64	25	46	53	27	20 (8")	7.6kg	130

Measurements in cm, Weight in kg