

# Fact Sheet

## Bathroom

### Shower Stool - Tool-less Assembly

The RM Compact Aluminium Shower Stool is a portable easy to assemble shower stool providing secure safe option for those who require a seat, without backrest, in the shower. A very simple click together frame; seat already attached, fixed armrests, ideal for permanent or temporary solution in the bathroom.












- \* Lightweight rust proof aluminium frame
- \* Durable powder coat finish
- \* Very easy to assemble
- \* Dismantle for storage / transportation
- \* Height adjustable
- \* Moulded plastic seat with drain holes
- \* Non-slip rubber tips
- \* Weight capacity 110kg



Product Code: RM392A-KD

#### Specifications

overall height	overall depth	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
								
66-78.5	47	54	22	47	47-59.5	19	2.15	110

Measurements in cm, Weight in kg