

Bathroom

Aluminium Shower Stool RM403










The RM lightweight aluminium shower stool with durable “S” bend side frame offering maximum comfort and space when seated. With large rubber suction ferrule to prevent slipping and provide maximum safety.



Product Code: RM403

- * Aluminium frame
- * Moulded plastic seat with drainage holes
- * “S” bend side frame provides maximum space
- * Large rubber suction ferrule to prevent slipping
- * Weight capacity 160kg

Specifications

overall height	overall depth	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
								
67-82	44	58	33	50	49-64	21	2.8	160

Measurements in cm, Weight in kg