

Bathroom

Solo Wall-mounted Seats

Solo wall-mounted stools and chairs provide a compact solution, offering a secure showering seat that can be folded away against the wall when not in use, ideal for areas where space is limited.



Product Code: VB544

- * All alloy lightweight frame, robust construction
- * Anti-bacterial finish to the ergonomic contoured padded PU seat
- * Fine adjustment leg for extra stability

Specifications

Seat dimensions36x28
 Seat height 47-54
 Weight capacity 159

- * Lightweight aluminium
- * Simple to install and maintain
- * Clip on, clip off seat with drainage for easy cleaning



Product Code: VB634

Specifications

Seat dimensions 49x38
 Seat height -stool46-54
 Seat height -chair48-61
 Width between arms.....46
 Depth from wall49
 Weight capacity..... 190

Measurements in cm, weight in kg



Product Code: V535