

## Bathroom



# Multi Purpose Stool with Rotating Seat

This versatile, multi-purpose stool can be used in a whole host of different settings around the home. With a comfortable rotating seat top making it ideal for those who need assistance when turning from side to side when seated.



Product Code: VB511R

- \* Durable, lightweight rust proof aluminium frame
- \* Suitable for use in bath / bathroom
- \* Rotating seat top
- \* Height Adjustable, 8 settings
- \* Suitable for children or adults
- \* Non-slip rubber tips
- \* Retail packaging
- \* Weight capacity 135kg

### Specifications

|                       |         |
|-----------------------|---------|
| Height (mm) .....     | 350-525 |
| Seat Depth (mm) ..... | 330     |
| Seat Width (mm) ..... | 330     |

