

## Wheelchair

### Strongback Excursion 12

The Strongback Excursion wheelchair is the first off-the-shelf standard folding wheelchair to provide frame-integrated effective posture support. Its unique ergonomic shape promotes healthy upright sitting and ensures comfort from the first moment. Unlike standard folding manual wheelchairs, the STRONGBACK Excursion conforms to and supports the natural and healthy lumbar curve of the spine thus preventing the user from slouching forward and squeezing spinal discs with a tremendous amount of unhealthy pressure - decreasing this pressure results in healthier discs, a stronger back, and long-lasting comfort while seated.

- \* Lightweight compact foldable design
- \* Easily adjustable posture support backrest
- \* Folding backrest
- \* Attendant brake
- \* Swing away detachable leg rests
- \* PU 12" rear wheel
- \* Weight capacity 135kg



**Product Code:** RMSB1001 40cm (16")  
 RMSB1012 50cm (20")

#### Specifications

overall width	seat width	armrest height	seat height	backrest height	overall height	seat depth	overall depth	unit weight	weight capacity
61 72	40cm 50cm	67 67	51 51	40 40	92 92	45 45	105 105	10 10	135 135