

Wheelchair

Strongback 24

The Strongback 24 is the first off-the-shelf standard folding manual wheelchair to provide frame-integrated effective posture support. Its unique ergonomic shape promotes healthy upright sitting and ensures comfort from the first moment. The Strongback 24 is as convenient, easy, and simple to operate as conventional wheelchairs. In addition to the outstanding ergonomics, the Strongback 24 also scores with stylishly curved lines and its elegant yet sporty color scheme.

- * Lightweight compact foldable design
- * Easily adjustable posture support backrest
- * Folding backrest
- * Attendant brake
- * Swing away detachable leg rests
- * Quick release 24" wheel
- * Weight capacity 135kg



Product Code: RMSB1006 40cm (16")
 RMSB1007 46cm (18")
 RMSB1010 50cm (20")

Specifications

| overall width | seat width | armrest height | seat height | backrest height | overall height | seat depth | overall depth | unit weight | weight capacity |
|---------------|------------|----------------|-------------|-----------------|----------------|------------|---------------|-------------|-----------------|
| 64 | 40cm | 67 | 50 | 40 | 94.5 | 45 | 102 | 12.9net11.2 | 135 |
| 69 | 45cm | 67 | 50 | 40 | 94.5 | 45 | 102 | 13 net11.4 | 135 |
| 74 | 50cm | 67 | 50 | 40 | 94.5 | 45 | 102 | 13.1net11.5 | 135 |