

Bathroom



Strood Aluminium Shower Stool

The Strood is a high quality shower stool designed for heavy duty use, with a maximum user weight of 159kg. Finished in durable hammertone powder coat helps prevent against scratches and scuffs.



Product Code: VB635AL

- * Lightweight rust proof aluminium frame
- * Durable hammertone powder coat finish
- * Height adjustable
- * Moulded plastic seat with drain holes
- * Clip off seat for easy clean
- * Non-slip rubber tips
- * Weight capacity 159kg

• **Australian Standards**
AS/NZS 3973:2009

Specifications

overall height	overall depth	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
56-72	42	52	28	41	46-62	15	2.8	159

Measurements in cm, Weight in kg