

Bathroom

Tri Corner Stool with Padded Seat

The Tri Corner Shower Stool is an excellent option when space is at a minimum. Lightweight, sturdy with rubber PVC stoppers for anti-slip stability. Contoured padded seat makes this chair a comfortable stool for use in the bathroom.



- * Ideal where space is limited
- * Lightweight aluminium (3.16kg)
- * Height adjustable
- * Drain holes under padded seat
- * Non-slip rubber tips
- * Weight capacity 136kg

Specifications

Height (mm)	490-650
Seat Depth (mm)	450
Seat Width (mm)	420
Maximum user weight (kg)	136

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com