

Wheelchair

Urban Wheelchair

Quality Foldable Transit Wheelchair makes life a bit simpler, for both the occupant and the attendant. A smart, sporty and classy looking compact unit, it is comfortable and user friendly.











- * Ultra light frame
- * Folding backrest
- * Attendant Brake
- * Swing up leg rests
- * Solid 6" front castor, 20" rear wheel
- * Rear wheel locking brake
- * Weight capacity 120kg



Product Code: WCT20

Specifications

Product Code: WCT16

overall width	seat width	armrest height	seat height	backrest height	overall height	seat depth	overall depth	unit weight	weight capacity
 16" - 58 20" - 64	 46 46	 23 23	 48 48	 41 41	 86 86	 40 40	 100 112	 10.5 11.5	 120 120

Measurements in cm, Weight in kg