

Pressure Care

Wheelchair Cushions

Foam / Memory Foam / Vinyl Cover



Specifications

Product Code: VA123SA 46x40x5cm
VA124SA 40x40x8cm
VA125SA 46x40x8cm

- * Quality foam with generous memory foam top layer
- * Super-soft vinyl cover with zip closure is easy to clean
- * Integrated ties secure the cushion to the wheelchair frame
- * Ideal for use in the car or around the home
- * Net weight 0.4kg

Foam / Memory Foam / Fabric Cover

- * Deluxe pressure relief orthopaedic wheelchair cushion
- * Convoluted (eggcrate) memory foam helps stimulate circulation
- * Allows air circulation, keeping you cool
- * 3 cm thick high density foam
- * 4.5 cm memory convoluted foam topper.
- * Integrated tie handles
- * Removable washable cover
- * Ideal for use in cars and around the home
- * Weight capacity 100kg
- * Net weight 0.75kg

Specifications

Product Code: VA126WA 43x43x8.5cm

