

Fact Sheet

Bathroom



Grab Rails

Grab rails are installed to provide greater independence, providing stability when using bathroom facilities and for general use around the home.

Solo Contract Hinged Arm Support

Product Code: VY423



- * Ideal for use in the bathroom where safety is a concern
- * Simple fold up / fold down mechanism with arm support leg
- * Epoxy coated ensuring durability
- * 750mm long bar

Solo Easigrip

450mm (18") VY394

* Durable epoxy coated steel

* Rounded ends for safety



Product Code: 300mm (12") VY393



Depth to wall (mm)	/60
Width (mm)	100
Leg Height (mm)	710-870
Leg from wall (mm)	555
Weight Capacity (kg)	115

President Ribbed Grab Bar

- * Durable UPVC finished in white
- * Ribbed grip for improved safety
- * Reinforced fixing points with cover caps
- * High quality AVS flanges
- * 32mm tube
- * Depth from wall 60mm



Product Code: 300mm (12") VY451

450mm (18") VY452



Product Code: 450mm (18") VY426

Product Code: 300mm (12") VY425