

## Bathroom



### Medina Bath Board

The Medina Bath Board is a high quality, reinforced plastic bath board.

- \* Modern design
- \* Slatted surface provides adequate drainage
- \* Wipe clean surface
- \* Simple to install with self locking cleats
- \* Both available with or without handle
- \* Weight Capacity 159kg



Product Code: VR116A (L:600mm)  
VR117A3 (L: 685mm)  
VR118A (L:710mm)

#### Specifications

- Overall Depth - 4 slat (mm) 250
- Overall Depth—6 slat (mm) 350
- Seat Height (mm) 140
- Seat Width (mm) 685
- Width Adjustment 535-625
- Net weight 3.2kg



Product Code: VR117H



Product Code: VR125A (L:685mm)

### Medina Bath Seat—Bariatric

The Medina Bath Seat provides extra height enabling easy access in and out of the bath.

- \* Constructed from high quality, fully reinforced slatted plastic
- \* Available in 2 heights (6" & 8")
- \* Suitable for majority of bath's
- \* Slatted for easy drainage
- \* Modern easy to clean surface
- \* Secured to bath with 4 large suction caps
- \* Easily removed by pulling on the integral tabs
- \* Finished in soft grey
- \* Weight Capacity 190kg

#### Specifications

- Dimensions (6") D290xW450xH150
- Dimensions (8") D290xW450xH200
- Net weight 2.4kg



Product Code: VR127A (6")  
VR128A (8")