

Bathroom



Raised Toilet Seats

Raised Toilet Seats are a convenient, ergonomically designed raised toilet seat. Designed to fit the majority of toilet bowl shapes. Manufactured from high quality ABS recyclable plastic. The side mounted brackets of the raised toilet seat have rubber over-moulds and the front skirt has a rubberised pad for maximum grip and stability. Also available padded 4".



Product Code: VR224C

Product Code: VR224D



Product Code: VR224H



Product Code: VR224F

- * Ergonomically designed for comfort
- * Smooth aperture interior for easy clean
- * Large aperture and front cutaway
- * 3 heights 2" (50mm) 4" (100mm) 6" (150mm)
- * Will fit most toilet suite
- * With or without lid



Product Code: VR224G



Product Code: VR223P



- * Padded for extra comfort
- * Removable cover can be washed
- * Secure Velcro fastening

Specifications

Aperture Width (mm).....225
Aperture Depth (mm).....270
ColourWhite
Seat Width (mm).....380
Seat Depth (mm).....400