

Bathroom



Stackable Bath Seat

The Plastic Bath Step features an extra-large platform with a textured, non-slip top surface and rubberised feet for increased safety. These handy bath steps can be securely stacked, raising the height in 2" increments, and allowing multiple units to be stacked and stored easily when not in use. The Plastic Bath Step is a versatile mobility aid that can also be used around the home, both inside and outside.



Product Code: VR278D

- * Versatile and lightweight
- * Secure, stackable step
- * Extra large platform
- * Non-slip textured surface
- * Rubber feet for safety
- * Weight capacity 120kg

Specifications

Height (mm)	100
Depth (mm)	420
Width (mm)	420
Net weight (kg)	1.05
Maximum user weight (kg)	120