

Beaumont Leg Rest (With Castors)

The top of the economical Beaumont Leg Rest has a flat surface, rather than a contoured one as featured on the Welling. The vinyl top of the Beaumont Leg Rest has heat sealed seams for improved hygiene and the frame is finished in a tough epoxy polymer finish. These padded leg rests are manufactured using a fire retardant, easy to clean fabric. Available with or without wheeled castors.



Features:

- Quality made leg rest ideal for use around the home to keep legs elevated comfortably
- Wipe clean vinyl finish
- Heat sealed seams
- Variable height (385 mm - 535 mm)

PRODUCT SPECIFICATION

Length (cm).....	310
Height (mm)	385-535
Width (mm)	500
Product Dimensions (mm)	385x500x310
Net weight (kg)	4.2
Packaged dimensions (HxWxL mm)	350x550x420
Display Footprint (mm)	385x500x310
Transport dimensions (mm)	350x550x420
Gross Weight (Packaged) (kg)	4.85
Packaging.....	Unpackaged
Item priced per	1
Outer box (HxWxL mm).....	550x420x330
Outer box quantity.....	2
Gross weight outer carton (kg).....	9.3
Colour.....	Black and brown
Configuration.....	With Castors
Height of Leg Rest (mm).....	385-535
Maximum User Weight (kg)	101
Maximum User Weight (stones).....	16