

E-TRAVELLER 140 CARBON

FOLDING ELECTRIC WHEELCHAIR



LITHIUM BATTERIES



FOLDS FOR EASY
TRANSPORT



LIGHTWEIGHT CARBON
FIBRE FRAME



E-TRAVELLER 140 CARBON

FOLDING ELECTRIC WHEELCHAIR

The E-Traveller 140 Carbon is a lightweight electric wheelchair made from carbon fibre that weighs just 17kg when the batteries are removed. An easy-folding mechanism makes this wheelchair ideal for frequent transport, either in the boot of your car, on a plane or on public transport. The Carbon is ideal for use in tighter confined areas as it has an unmatched turning circle.

SPECIFICATIONS

Weight Capacity	140kg	Battery Weight	1.45kg x 2
Dimensions LxWxH (mm)	950x590x940	Max Slope	8°
Folded Dimensions (mm)	590x275x760	Tyres	Solid - Puncture Proof
Max Speed	6kph	Drive Wheel	203mm
Battery Range	Up to 20km	Motor	2 x 180W
Battery	24V/6ah Lithium x 2	Joystick Control	Yes
Total Weight Without Batteries	17kg	Total Weight With Batteries	20.45kg



**SCAN OR CLICK
TO FIND OUT MORE**