




- Lightweight, easy to use skin protection for up to high risk and including grade 2 damage.
- Ultra soft, touch-and-close adjustable fastening material to avoid skin damage
- Dark grey colour to maintain a smart, clinical appearance
- Generous inspection hole to allow easy heel inspection
- Shaped foam interior helps to maintain airflow to reduce heat build up, maximise contact area and thus minimise interface pressures.
- Ideal for both heels and elbows.
- Neatly shrink wrapped in pairs for easy storage.
- Manufactured in Derbyshire. 

NB: Not Suitable for ambulation. Single patient use only

Product Code	Quantity
HB19030	30 Pairs (box)
HB19020	20 Pairs (box)
HB190	1 Pair