

KYMCO KOMFY 8

MEDIUM SIZED MOBILITY SCOOTER



**FULLY ADJUSTABLE
SUSPENSION**



**COMFORTABLE
PNEUMATIC TYRES**



**USB PHONE
CHARGING PORT**



KYMCO KOMFY 8

MEDIUM SIZED MOBILITY SCOOTER

The Kymco Komfy 8 is designed to be user friendly and features a dashboard with large controls, delta grip handlebars, adjustable mirrors and a handy USB phone charging port. A full LED lighting package, adjustable captain style seat, front and rear suspension tick all the boxes. Both the front coil over spring suspension and the rear mono-shock absorber are fully adjustable, combined with the 10-inch pneumatic tyres the Kymco Komfy 8 is packed with many great features for a medium sized scooter.

SPECIFICATIONS

Weight With Batteries	93kg	Tyres	Pneumatic
Battery Weight	-	Dimensions LxW (mm)	1280x600
Heaviest Part	-	Swivel Seat	Yes
Weight Capacity	160kg	Distance Between Armrests	460mm-660mm
Adjustable Tiller Angle	Yes	Drive Wheel	254mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	45mm
Battery	12V/35ah x 2	Seat Dimensions	460mm x 440mm
Battery Range	Up To 32km	Seat Height From Ground	580mm-640mm
Maximum Speed	10km/h	Slope Ability	10°
Motor	500W 4-Pole	Turning Circle	1560mm
Number of Components	1	Suspension	Adjustable Front and Rear

Colours Subject to Availability



Dealer Details