

# KYMCO MAXER

LARGE SIZED MOBILITY SCOOTER



**FRONT AND REAR  
LIGHTING**



**WATERPROOF  
BACK LIT LCD DISPLAY**



**FULLY ADJUSTABLE  
SUSPENSION**



# KYMCO MAXER

## MEDIUM SIZED MOBILITY SCOOTER

The Maxer is one of the largest scooters in our range and is made to tackle all conditions. Yet, its slim line design will allow you to easily navigate doorways and other entries, while the large 700W 4-pole motor offers enough power to navigate the more challenging terrains safely. The adjustable suspension makes it perfect for most outdoor use, and the 80Ah battery will allow it to travel up to 57km on a single charge. A fully adjustable, contoured captain style seat will provide the ultimate in comfort and versatility.

### SPECIFICATIONS

<b>Weight With Batteries</b>	138kg	<b>Tyres</b>	Pneumatic
<b>Battery Weight</b>	-	<b>Dimensions LxW (mm)</b>	1420x700
<b>Heaviest Part</b>	-	<b>Swivel Seat</b>	Yes
<b>Weight Capacity</b>	200kg	<b>Distance Between Armrests</b>	480mm-680mm
<b>Adjustable Tiller Angle</b>	Yes	<b>Drive Wheel</b>	335mm
<b>Anti-Tip Mechanism</b>	Rear Anti-Tip Wheels	<b>Ground Clearance</b>	100mm
<b>Battery</b>	12V/80ah x 2	<b>Seat Dimensions</b>	480mm x 460mm
<b>Battery Range</b>	Up To 57km	<b>Seat Height From Ground</b>	650mm-710mm
<b>Maximum Speed</b>	10km/h	<b>Slope Ability</b>	12°
<b>Motor</b>	700W 4-Pole	<b>Turning Circle</b>	2125mm
<b>Number of Components</b>	1	<b>Suspension</b>	Adjustable Front and Rear

Colours Subject to Availability



Dealer Details