

# QUINGO FLYTE

PORTABLE 5 WHEEL MOBILITY SCOOTER



**INNOVATIVE  
5-WHEEL DESIGN**



**SELF LOADING  
DOCKING STATION**



**EASY TO USE  
REMOTE**



# QUINGO FLYTE

## SELF LOADING 5 WHEELED SCOOTER

The Self-loading portable and perfect solution for those who struggle with lifting. Ideal indoors or out and loads into most hatchbacks/SUVs in under 60 seconds. The flyte just got even better with the lighter (the heaviest component weighs only 8.9kg), easier to fit & more flexible all-new MK2 docking station solution from most hatchbacks and SUVs with the SUV adaptor kit.

### SPECIFICATIONS

<b>Weight With Batteries</b>	116.5kg	<b>Tyres</b>	Pneumatic
<b>Battery Weight</b>	-	<b>Dimensions LxW (mm)</b>	1195x620
<b>Heaviest Part</b>	-	<b>Swivel Seat</b>	No
<b>Weight Capacity</b>	159kg	<b>Distance Between Armrests</b>	460mm-660mm
<b>Adjustable Tiller Angle</b>	Yes	<b>Drive Wheel</b>	250mm
<b>Anti-Tip Mechanism</b>	Rear Anti-Tip Wheels	<b>Ground Clearance</b>	75mm
<b>Battery</b>	12V 22Ah x 4	<b>Seat Dimensions</b>	410mm x 400mm
<b>Battery Range</b>	Up To 37km	<b>Seat Height From Ground</b>	710mm
<b>Maximum Speed</b>	6.4km/h	<b>Slope Ability</b>	6°
<b>Motor</b>	250W 2-Pole	<b>Turning Circle</b>	1020mm
<b>Number of Components</b>	1 + Docking Station	<b>Suspension</b>	No