

# QUINGO PLUS

5 WHEEL MOBILITY SCOOTER



**INNOVATIVE  
5-WHEEL DESIGN**



**LARGE  
COMFORTABLE SEAT**



**KERB CLIMBING  
TECHNOLOGY**



# QUINGO PLUS

## 5 WHEEL MOBILITY SCOOTER

The Quingo Plus is an upgraded model of the Quingo classic, featuring everything we love about the Classic and the 5-wheel Quingo design plus additional comfort and usability features, such as high and low-speed selection, a broader range and a large adjustable captain seat. The Plus is classed as a medium-sized mobility scooter but can outperform many other large gophers. However, thanks to the 5-wheel design, it still allows for the same maneuverability that you are used to from smaller 4-wheel models.

### SPECIFICATIONS

<b>Weight With Batteries</b>	103kg	<b>Tyres</b>	Pneumatic
<b>Battery Weight</b>	-	<b>Dimensions LxW (mm)</b>	1270x650
<b>Heaviest Part</b>	-	<b>Swivel Seat</b>	Yes
<b>Weight Capacity</b>	159kg	<b>Distance Between Armrests</b>	450mm-630mm
<b>Adjustable Tiller Angle</b>	Yes	<b>Drive Wheel</b>	250mm
<b>Anti-Tip Mechanism</b>	Rear Anti-Tip Wheels	<b>Ground Clearance</b>	55mm
<b>Battery</b>	12V 50Ah x2	<b>Seat Dimensions</b>	480mm x 440mm
<b>Battery Range</b>	Up To 48km	<b>Seat Height From Ground</b>	600mm - 700mm
<b>Maximum Speed</b>	6.4km/h	<b>Slope Ability</b>	12°
<b>Motor</b>	400W 4-Pole	<b>Turning Circle</b>	1070mm
<b>Number of Components</b>	1	<b>Suspension</b>	Yes