

# QUINGO TOURA 2

5 WHEEL MOBILITY SCOOTER



**INNOVATIVE  
5-WHEEL DESIGN**



**DIGITAL  
DASHBOARD**



**FRONT BASKET  
ATTACHMENT**



# QUINGO TOUREA 2

## 5 WHEEL MOBILITY SCOOTER

The Quingo Tourea 2 is the top of the line scooter of the 5-wheel Quingo range. It features countless comfort and safety features, all while maintaining excellent maneuverability for such a substantial machine. With a carrying capacity of 220kg and a travel range of up to 88km this scooter makes the perfect car replacement. Featuring an easy to read digital dash with a battery meter gauge, easy to use speed dial and push-button lighting control the dash has everything that you need exactly where you need it.

### SPECIFICATIONS

<b>Weight With Batteries</b>	163.9kg	<b>Tyres</b>	Pneumatic
<b>Battery Weight</b>	-	<b>Dimensions LxW (mm)</b>	1420x700
<b>Heaviest Part</b>	-	<b>Swivel Seat</b>	Yes
<b>Weight Capacity</b>	220kg	<b>Distance Between Armrests</b>	540mm-720mm
<b>Adjustable Tiller Angle</b>	Yes	<b>Drive Wheel</b>	305mm
<b>Anti-Tip Mechanism</b>	No	<b>Ground Clearance</b>	100mm
<b>Battery</b>	12V 100Ah x2	<b>Seat Dimensions</b>	540mm x 470mm
<b>Battery Range</b>	Up To 88km	<b>Seat Height From Ground</b>	720mm - 790mm
<b>Maximum Speed</b>	10km/h	<b>Slope Ability</b>	12°
<b>Motor</b>	650W 4-Pole	<b>Turning Circle</b>	1320mm
<b>Number of Components</b>	1	<b>Suspension</b>	Yes