

# QUINGO VITESS 2

5 WHEEL MOBILITY SCOOTER



**INNOVATIVE  
5-WHEEL DESIGN**



**DIGITAL  
DASHBOARD**



**REAR ANTI-TIP  
WHEELS**



# QUINGO VITESS 2

## 5 WHEEL MOBILITY SCOOTER

The Quingo Vites 2 is the perfect mobility solution for more hilly terrain. It features a Digital Dash, superior lighting & more power. The unique 5 wheel stability system on all Quingo's provides the excellent posture and driving position of a 3 wheeler with the stability and safer cornering dynamics of a 4 wheeler. Along with the best possible anatomical driving position all Quingo's also benefit from a staggering 80% more foot area than in a comparable length four-wheel vehicle.

### SPECIFICATIONS

<b>Weight With Batteries</b>	116kg	<b>Tyres</b>	Pneumatic
<b>Battery Weight</b>	-	<b>Dimensions LxW (mm)</b>	1420x700
<b>Heaviest Part</b>	-	<b>Swivel Seat</b>	Yes
<b>Weight Capacity</b>	182kg	<b>Distance Between Armrests</b>	450mm-630mm
<b>Adjustable Tiller Angle</b>	Yes	<b>Drive Wheel</b>	178mm
<b>Anti-Tip Mechanism</b>	Rear Anti-Tip Wheels	<b>Ground Clearance</b>	75mm
<b>Battery</b>	12V 50Ah x2	<b>Seat Dimensions</b>	470mm x 430mm
<b>Battery Range</b>	Up To 48km	<b>Seat Height From Ground</b>	640mm - 740mm
<b>Maximum Speed</b>	10km/h	<b>Slope Ability</b>	12°
<b>Motor</b>	400W 4-Pole	<b>Turning Circle</b>	1130mm
<b>Number of Components</b>	1	<b>Suspension</b>	Yes