



	<b>The Select Range</b> - In Studio Vinyls Available in Single Motor	<b>The Signature Range</b> - In Angora Fabrics Available in Spacesaver Single Motor & Infinite Position Dual Motor		<b>The Luxury Range</b> - Genuine Hide Leather Available in Spacesaver Single Motor & Infinite Position Dual Motor		<b>The Chadwick Range</b> - Soft Touch Fabric Available in Quad Motors	<b>The Chadwick Range</b> - Leather Available in Quad Motors	<b>The Chadwick Range</b> - Oxford Plush Leather Available in Quad Motors
Motor	Select - Single Motor	Signature - Single Motor	Signature - Dual Motor	Luxury - Single Motor	Luxury - Dual Motor	Infinite - Quad Motor	Infinite - Quad Motor	Infinite - Quad Motor
Overall Height	105cm	105cm	105cm	105cm	105cm	110cm	110cm	110cm
Overall Width	75cm	75cm	75cm	75cm	75cm	85cm	85cm	85cm
Overall Depth	75cm	75cm	75cm	75cm	75cm	80cm	80cm	80cm
Armrest Height (From floor)	62cm	62cm	62cm	62cm	62cm	64cm	64cm	64cm
Seat Height (From floor)	48cm	48cm	48cm	48cm	48cm	50cm	50cm	50cm
Seat Width (Measured on sitting area)	52cm	52cm	52cm	52cm	52cm	56cm	56cm	56cm
Seat Depth (Measured on sitting area)	50cm	50cm	50cm	50cm	50cm	50cm	50cm	50cm
Seat Max. Height (Lifted)	56cm	58cm	56cm	58cm	57cm	57cm	57cm	57cm
Backrest Length (Fully reclined)	32cm	32cm	60cm	32cm	60cm	60cm	60cm	60cm
Max. Legrest Height (Lifted)	46cm	46cm	48cm	46cm	48cm	52cm	52cm	52cm
Safe Working Load (SWL)	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs
Foam Density	28/70	28/70	28/70	28/70	28/70	35/70	35/70	35/70
Upholstery	Studio Vinyl	Premium Angora Fabrics	Premium Angora Fabrics	Genuine Hide Leathers Bovine Yearling	Genuine Hide Leathers Bovine Yearling	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)
Colour Options	Charcoal black	Symphony Blue	Symphony Blue	Charcoal Black	Charcoal Black	Indigo Blue	Taupe Beige	Colonial Hemp Beige
	Mocha Beige	Dark Coffee Bean	Dark Coffee Bean	Espresso Brown	Espresso Brown	Oyster Latte	Espresso Dark Brown	Jet Carbon Black
	Mushroom Brown	Mushroom Grey	Mushroom Grey	Taupe Beige	Taupe Beige	Coffee Bean Brown		Dark Cigar Brown
	Chocolate Brown							Marine Blue

The Chadwick and Mayfair Lift Chair collections are ILS's premium lift chair range. They come in wide range of colours and materials to suit everyone's budget.

**"BEST SEAT IN THE HOUSE"**