















































	The Select Range - In Studio Vinyls Available in Single Motor	The Signature Range - In Angora Fabrics Available in Spacesaver Single Motor & Infinite Position Dual Motor		The Luxury Range - Genuine Hide Leather Available in Spacesaver Single Motor & Infinite Position Dual Motor		The Chadwick Range - Soft Touch Fabric Available in Quad Motors	The Chadwick Range - Leather Available in Quad Motors	The Chadwick Range - Oxford Plush Leather Available in Quad Motors
								
Motor	Select - Single Motor	Signature - Single Motor	Signature - Dual Motor	Luxury - Single Motor	Luxury - Dual Motor	Infinite - Quad Motor	Infinite - Quad Motor	Infinite - Quad Motor
 Overall Height	105cm	105cm	105cm	105cm	105cm	110cm	110cm	110cm
 Overall Width	75cm	75cm	75cm	75cm	75cm	85cm	85cm	85cm
 Overall Depth	75cm	75cm	75cm	75cm	75cm	80cm	80cm	80cm
 Armrest Height (From floor)	62cm	62cm	62cm	62cm	62cm	64cm	64cm	64cm
 Seat Height (From floor)	48cm	48cm	48cm	48cm	48cm	50cm	50cm	50cm
 Seat Width (Measured on sitting area)	52cm	52cm	52cm	52cm	52cm	56cm	56cm	56cm
 Seat Depth (Measured on sitting area)	50cm	50cm	50cm	50cm	50cm	50cm	50cm	50cm
 Seat Max. Height (Lifted)	56cm	58cm	56cm	58cm	57cm	57cm	57cm	57cm
 Backrest Length (Fully reclined)	32cm	32cm	60cm	32cm	60cm	60cm	60cm	60cm
 Max. Legrest Height (Lifted)	46cm	46cm	48cm	46cm	48cm	52cm	52cm	52cm
 Safe Working Load (SWL)	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs
Foam Density	28/70	28/70	28/70	28/70	28/70	35/70	35/70	35/70
Upholstery	Studio Vinyl	Premium Angora Fabrics	Premium Angora Fabrics	Genuine Hide Leathers Bovine Yearling	Genuine Hide Leathers Bovine Yearling	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)
Colour Options	 Charcoal black	 Symphony Blue	 Symphony Blue	 Charcoal Black	 Charcoal Black	 Indigo Blue	 Taupe Beige	 Colonial Hemp Beige
	 Mocha Beige	 Dark Coffee Bean	 Dark Coffee Bean	 Espresso Brown	 Espresso Brown	 Oyster Latte	 Espresso Dark Brown	 Jet Carbon Black
	 Mushroom Brown	 Mushroom Grey	 Mushroom Grey	 Taupe Beige	 Taupe Beige	 Coffee Bean Brown		 Dark Cigar Brown
	 Chocolate Brown							 Marine Blue

The Chadwick and Mayfair Lift Chair collections are ILS's premium lift chair range. They come in wide range of colours and materials to suit everyone's budget.

"BEST SEAT IN THE HOUSE"

