

SOLAX MALETA V2

AUTOMATIC FOLDING MOBILITY SCOOTER



DIGITAL DASHBOARD



**REMOTE CONTROL
AUTOMATIC FOLDING**



**LIGHTWEIGHT
COMPACT DESIGN**



SOLAX MALETA V2

AUTOMATIC FOLDING MOBILITY SCOOTER

The Maleta V2 is packed with new innovative modern design features. The V2 sees the addition of a charging port on the dashboard, a remote control for folding, headlights, ambient deck lights and ergonomic horn shaped handle bars. It is safe, comfortable, easy to operate, convenient to be folded and carried, lightweight and steady. The frame is made of a lightweight aluminium alloy. It can be folded and unfolded without any tools or effort. It has 2 automatic folding functions and a manual safety fold override. A lightly padded seat and spring style backrest will provide you with support during journeys.

SPECIFICATIONS

Weight With Batteries	25.55kg	Tyres	Solid - Puncture Proof
Battery Weight	2 x 0.85kg	Dimensions LxWxH (mm)	970x460x950
Weight Without Batteries	23.85kg	Folded Dimensions (mm)	280x460x650
Weight Capacity	125kg	Lights	Front, Rear and Ambient Deck
Adjustable Tiller Height	Yes	Drive Wheel	178mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	40mm
Battery	24V/5.6ah Lithium x 2	Seat Dimensions	400mm x 315mm
Battery Range	Up To 15km	Seat Height From Ground	560mm
Maximum Speed	6km/h	Slope Ability	10°
Motor	120W	Turning Circle	1200mm
Number of Components	1 + Batteries	Suspension	No