

## SPECIFICATIONS

	TRANSPORT	STANDARD	HEAVY DUTY/ EXTRA WIDE	PLUS™ EXTRA WIDE	LIGHTWEIGHT	HIGH STRENGTH LIGHTWEIGHT
<b>SEAT WIDTH</b>	17, 19 in.	16 & 18 in.	20, 22 & 24 in.	26, 28 & 30 in.	16, 18 & 20 in.	16, 18 & 20 in.
<b>SEAT DEPTH</b>	16 inches	16 inches	18 inches 16 inches*	20 inches	16 inches	16 inches
<b>OVERALL WIDTH</b> FIXED ARM REMOVABLE ARM	21 1/2, 23 1/2 in. N/A	22 1/2, 24 1/2 in. 24 1/2, 26 1/2 in.	N/A 29 1/2, 30 1/2, 32 1/2 in.	N/A 35, 37, 39 in.	N/A 23, 25, 27 in.	N/A 23, 25, 27 in.
<b>OVERALL DEPTH</b> (W/O RIGGINGS)	24 inches	Fixed Arm: 32 in. Removable Arm 29 in.	32 inches	37 inches	30 1/2 in.	30 1/2 in.
<b>SEAT TO FLOOR HEIGHT</b> ADULT HEMI	19 inches N/A	19 3/4 in. 17 3/4 in.	19 3/4 in. 17 3/4 in.*	19 3/4 in. N/A	19 3/4 in. 17 3/4 in.	19 3/4 in. 17 3/4 in.
<b>BACK STYLE</b>	Fixed & Fold Down	Fixed	Fixed	Fixed	Fixed	Adjustable
<b>BACK HEIGHT</b>	17 inches	16 inches	16 inches	18 in.	16 inches	16 - 18 in.
<b>ARM STYLES</b>	Fixed Full Length	Fixed or Removable Full or Desk Length	Removable Full Length Removable Desk Length	Removable Desk Length	Removable Desk Length	Flip Back Removable Desk Length
<b>FRONT RIGGING</b>	Swing-Away Footrest	Fixed Footrest Swing-Away Footrest Elevating Legrest	Swing-Away Footrest Elevating Legrest	Swing-Away Footrest Elevating Legrest	Swing-Away Footrest Elevating Legrest	Swing-Away Footrest Elevating Legrest
<b>REAR AXLE</b>	Single Position	Dual Position	Single Position Dual Position*	Single Position	Dual Position	Dual Position
<b>REAR WHEELS</b>	8x1 in. Rubber	Composite Mag Style	Composite Mag Style	Steel Mag Style	Composite Mag Style	Composite Mag Style
<b>HAND RIMS</b>	N/A	Composite	Chromed Steel	Chromed Steel	Composite	Aluminum
<b>WHEEL LOCKS</b>	Push-to-Lock	Push-to-Lock	Push-to-Lock	Pull-to-Lock	Push-to-Lock	Pull-to-Lock
<b>FRONT CASTERS</b>	8x1 in. Rubber	8x1 in. Rubber	8x2 in. Rubber	8x2 in. Rubber	8x1 in. Rubber	8x1 in. Rubber
<b>UPHOLSTERY</b>	Nylon	Vinyl	Vinyl	Nylon	Nylon	Nylon

\*20 inch only