

Bathroom *Bariatric*

Bath Transfer Bench

The durable aluminium Bath Transfer Bench is placed over the bath for those who have a shower over the bath situation. The wider seat allows the client to slide across from the outside of the bath. Two large suction feet for secure placement. Fully height adjustable will fit over most baths.



- * Durable aluminium frame
- * Detachable back rest, can be attached to either side
- * Heavy duty suction tips
- * Detachable arm rest
- * Weight capacity 200kg

- **Complies with Australian Standards
AS/NZS 3973:2009**

Specifications

Overall height..... 84-99
Overall depth.....55
Overall width.....78
Seat depth41
Seat width68

Tub clearance..... 41-55
Armrest height30
Unit weight.....6
Weight capacity.....200

Measurements in cm, weight in kg

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com