

Lightweight Aluminium Wheelchair

The Lightweight Aluminium Wheelchair is available in both self propel and transit options. With many of the features of other chairs the lightweight wheelchair offers ease of movement for the user and easy folding for trouble-free transportation.



- * Lightweight Transit
- * Folding backrest (for storage or transport)
- * Care Assist brakes
- * Flip up armrest for side transfer
- * 12" rear wheel, 6" front castor
- * Swing away and removable leg rests
- * Solid tyres (no maintenance required)
- * Only available in red
- * Weight capacity 110kg

Specifications

Product Code: WH912-46

	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Footrest Height	Unit Weight	Weight Capacity
WH912-46	70	106	63	40	46	50	20	43	33-44	14.3	110

Measurements in cm, Weight in kg

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com