

Peak Over Bed Table

The Over Bed Table is a versatile mobile option for the around the house especially great for use in bedrooms. With a height adjustable tilting table top this is great for both reading and writing.

- * Height adjustable
- * 4 lockable 2" castors
- * Easy clean
- * Tilting top for reading and writing
- * Steel frame
- * Table top weight capacity 15kg



Specifications

Overall height.....	75-111	Under bed clearance.....	9
Overall depth	39	Unit weight.....	7.2
Overall width.....	68	Weight capacity	15
Top depth.....	39		
Top width	64		

Measurements in cm, weight in kg

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com