

Bathroom



Pembury Corner Stool

The Pembury height adjustable corner shower stool, an excellent option when space is at a premium. Featuring fine leg adjustment maximises stability and safety when using the shower. A clip on/clip off seat with drainage holes allows for easy cleaning. Constructed from lightweight aluminium makes it easy to manage.



- * Ideal where space is at a premium
- * Powder coated
- * Height Adjustable
- * Non-slip rubber tips
- * Weight capacity 159kg

Specifications

Height (mm)	470-620
Seat Depth (mm)	265
Seat Width (mm)	530
Overall Depth.....	380
Maximum user weight (kg)	159