

Ramps

Lightweight Suitcase

- * Suits majority of scooters and wheelchairs
- * Provides access to vehicles and buildings
- * Safety all-weather tread
- * Easy to use and transport

Specifications

Folded width 370mm
 Channel width..... 700mm
 Weight capacity 272kg (WCR180).....300kg



Product Code:

VA142C 910mm(3ft) net weight 9kg
VA142K 1520mm(5ft)net weight 14kg
WCR180 1830mm(6ft) net weight 19kg

Ramp Roll Up



- * Easy access to houses and cars
- * Portable & Compact
- * Supplied with handy carry bag (20x20x90cm)

Specifications

Length 910mm
 Width760mm
 Weight capacity300kg
 Net weight..... 9kg

Product Code: (3ft) **VA141**

Telescopic Channel

- * Gritted surface for increased grip
- * Lightweight

Specifications

Length collapsed (VA147F) 744mm
 (VA147H)835mm
 (VA147M)985mm
 (VA147S)1201mm
 Channel width189mm
 Weight capacity180kg



Product Code:

VA147F 1220mm(4ft) net weight 6kg
VA147H 1520mm(5ft) net weight 7kg
VA147M 1830mm(6ft) net weight 8kg
VA147S 2130mm(7ft) net weight 10kg



Folding Suitcase

- * Non slip surface
- * Tri-fold for easy storage and transport

Specifications

Folded width360mm
 Folded length 630 &915mm
 Channel width740mm
 Weight capacity 272kg

Product Code:

VA143D 1220mm(4ft) net weight 11kg
VA143S 1830mm(6ft) net weight 20kg