

Bathroom

Toilet Surrounds

Height adjustable Toilet Surrounds assist in getting on and off the toilet without having to install wall rails. Durable aluminium powder coated finish.



- * Height adjustable
- * Width adjustable (VR202)
- * Double transform bar for added strength
- * Coated with a tough white polymer anti-scratch finish
- * Moulded support arms provide improved grip
- * Will fit most toilet suites
- * Weight Capacity 190kg (VR205 & VR202)



- * Aluminium frame
- * Height adjustable hand grip
- * Detachable frame
- * Rubber tips to reduce slip
- * Weight Capacity 130kg



Specifications VR202 & VR205

Depth (mm)	450
Height (mm).....	660
VR202 int width (mm)	490-630
VR205 int width (mm)	484
Width (mm).....	580
Net weight (kg).....	4.29
Colour.....	White
Weight Capacity	190kg

Specifications RM418

Length (mm).....	520
Height (mm)	650-771
Width (mm).....	580
Net weight (kg)	2.6
Colour.	White
Weight Capacity	130kg