

Bathroom

Toilet Surround

Height adjustable Toilet Surrounds assist in getting on and off the toilet without having to install wall rails. Durable aluminium powder coat finish.

- * Aluminium frame
- * Height adjustable hand grip
- * Detachable frame
- * Rubber tips to reduce slip
- * Available in 54cm & 60cm width



Specifications RM418

Overall height	65-77
Overall depth	52
Overall width	58/62
Unit weight	2.6
Weight capacity	130

Measurements in cm, weight in kg

Ausnew Home Care Service

Suite 1/88 Rooty Hill Rd North, Rooty Hill NSW 2766 Australia

ABN : 31 620 493 941

Phone: 291596183

Email support@ausnewhomecare.com